


# COG




- 3D Multi directional COG
- One lined
- Specialized in Lifting
- For all parts of face and body

Thread type	Needle Size Gauge	Needle Length mm	Thread Length mm	Thread Thickness USP	Package Q.ty Pcs/pouch
Bidirectional 3D Cog (Cut Type, Blunt Needle)	19G	40mm	80mm	0-0	20
	19G	60mm	100mm	0-0	20
	19G*	100mm	150mm	0-0	20
Bidirectional 3D Cog (Cut Type)	21G	40mm	80mm	2-0	20
	21G	50mm	90mm	2-0	20
	21G*	60mm	100mm	2-0	20
	21G	90mm	150mm	2-0	20
	23G	40mm	80mm	3-0	20
	23G	50mm	90mm	3-0	20
	23G	60mm	100mm	3-0	20
Bidirectional 3D Cog (Non-Cut Type)	21G	40mm	50mm	2-0	20
	21G	50mm	70mm	2-0	20
	21G	60mm	90mm	2-0	20
	21G	90mm	150mm	2-0	20
	23G	40mm	50mm	3-0	20
	23G	50mm	70mm	3-0	20
	23G*	60mm	90mm	3-0	20
	23G*	90mm	150mm	3-0	20


\*Always available

# JALOR THREADS "Sweet Arrow"

## EXPECTED EFFICACY



- Non-heat-treated manufacturing method
- Maintains original strength of thread
- Long lasting effect of 1~2 years
- High anchoring strength
- High tensile strength
- Minimal tissue reaction
- Minimal trauma to tissue
- Wide area of collagen formation




## SPECIFICATIONS

Needle Type	Needle Diameter	Needle Length	Total Length	Barb Distance
L-Cannula type	18G	100mm	185mm	90mm
L-Cannula type	19G	100mm	185mm	90mm
L-Cannula type	18G	70mm	162mm	67mm
L-Cannula type	19G	70mm	162mm	67mm



# JALOR THREADS



# LIFTING USING PDO THREAD IMPLANTS

“Jalor threads” is a new lifting method based on the implanting of multiple Polydioxanone (PDO) mini threads, creating a vectorial network that acts as support for the tissue providing an immediate visible lifting effect while also favouring natural collagen production.

Treatment is relatively simple and fast (15-30 minutes). It involves the insertion of small needles loaded with Polydioxanone (PDO) mini threads that remain in the skin once the needles are removed.

The effectiveness of the treatment is clearly visible after 2-3 weeks, in fact day by day the patient will be able to see on its skin a significant improvement in the complexion and skin texture. The PDO threads are gradually reabsorbed over a period of 6 to 8 months via biological absorption or hydrolysis but the results can last between one and two years, depending on the characteristics of each patient. In fact the effect is not only connected to the time of permanence but also to the endogenous stimulation promoted by the threads. They produce an increase in collagen around the threads. Although the effects are immediately apparent, maximum tension is reached some time later (between two to three months).


It is advisable to carry out maintaining treatments in order to provide continuity to the results. The treatment stimulates skin rejuvenation, bringing vitality, elasticity and brighter skin tone. Treatment can be carried out at the Doctor's surgery with only local anaesthetic required. No hospitalisation required. No allergic reactions or rejection. No side effects. The availability of a variety of different needle sizes facilitates treatment and adaptation to distinct areas and patient characteristics.

It also has therapeutic effects, in fact improves tissue and enhancing blood and lymph circulation, reduce skin stress and slows down aging, delaying the biological clock, regenerate naturally the production of collagen and elastin, proving a more rested appearance.

**Excellent results in both face and body lifts.**




## MONO



- One lined basic type thread
- For lifting and volume up
- For all parts of face and body

Thread type	Needle Size Gauge	Needle Length mm	Thread Length mm	Thread Thickness USP	Package Q.ty Pcs/pouch
Mono	25G	40mm	50mm	5-0	50
Mono	25G	50mm	70mm	5-0	50
Mono	25G	60mm	90mm	5-0	50
Mono	25G	90mm	150mm	5-0	50
Mono	26G	40mm	50mm	5-0	50
Mono	26G	50mm	70mm	5-0	50
Mono	26G	60mm	90mm	5-0	50
Mono*	26G	90mm	150mm	5-0	50
Mono	27G	40mm	50mm	5-0	50
Mono	27G	50mm	70mm	5-0	50
Mono	27G	60mm	90mm	5-0	50
Mono	29G	25mm	30mm	6-0	50
Mono*	29G	40mm	50mm	6-0	50
Mono	29G	50mm	70mm	6-0	50
Mono*	29G	60mm	90mm	6-0	50
Mono*	30G	25mm	30mm	7-0	50
Mono	30G	40mm	50mm	7-0	50
Mono	30G	50mm	70mm	7-0	50


## SCREW



- One lined thread
- Screwed on the needle
- For lifting and volume up
- For all parts of face and body

Thread type	Needle Size Gauge	Needle Length mm	Thread Length mm	Thread Thickness USP	Package Q.ty Pcs/pouch
Screw	25G	40mm	50mm	5-0	50
Screw	25G	50mm	70mm	5-0	50
Screw	25G	60mm	90mm	5-0	50
Screw*	25G	90mm	150mm	5-0	50
Screw*	26G	40mm	50mm	5-0	50
Screw	26G	50mm	70mm	5-0	50
Screw	26G	60mm	90mm	5-0	50
Screw	26G	90mm	150mm	5-0	50
Screw	27G	40mm	50mm	5-0	50
Screw	27G	50mm	70mm	5-0	50
Screw	27G	60mm	90mm	5-0	50
Screw	29G	25mm	30mm	6-0	50
Screw	29G	40mm	50mm	6-0	50
Screw*	29G	50mm	70mm	6-0	50
Screw	29G	60mm	90mm	6-0	50
Screw	30G	25mm	30mm	7-0	50
Screw	30G	40mm	50mm	7-0	50
Screw	30G	50mm	70mm	7-0	50


## TWIN



- Two lined and twisted thread
- For lifting and volume up
- For all parts of face and body

Thread type	Needle Size Gauge	Needle Length mm	Thread Length mm	Thread Thickness USP	Package Q.ty Pcs/pouch
Twin	25G	40mm	50mm	6-0	50
Twin	25G	50mm	70mm	6-0	50
Twin	25G	60mm	90mm	6-0	50
Twin	25G	90mm	150mm	6-0	50
Twin	26G	40mm	50mm	7-0	50
Twin*	26G	50mm	70mm	7-0	50
Twin	26G	60mm	90mm	7-0	50
Twin*	26G	90mm	150mm	7-0	50
Twin	27G	40mm	50mm	7-0	50
Twin	27G	50mm	70mm	7-0	50
Twin	27G	60mm	90mm	7-0	50

## TWIN SCREW



- Two lined and twisted thread
- Screwed on the needle
- For lifting and volume up
- For all parts of face and body

Thread type	Needle Size Gauge	Needle Length mm	Thread Length mm	Thread Thickness USP	Package Q.ty Pcs/pouch
Twin Screw	25G	40mm	50mm	6-0	50
Twin Screw	25G	50mm	70mm	6-0	50
Twin Screw	25G	60mm	90mm	6-0	50
Twin Screw	25G	90mm	150mm	6-0	50
Twin Screw	26G	40mm	50mm	7-0	50
Twin Screw*	26G	50mm	70mm	7-0	50
Twin Screw	26G	60mm	90mm	7-0	50
Twin Screw	26G	90mm	150mm	7-0	50
Twin Screw	27G	40mm	50mm	7-0	50
Twin Screw	27G	50mm	70mm	7-0	50
Twin Screw*	27G	60mm	90mm	7-0	50

\*Always available